


# DROP ZONE™ the CLOT INSIDE

## DROP ZONE TECHNOLOGY

taking clot capture to the next level



Nevala™

Designed for **1<sup>st</sup> PASS SUCCESS** with **ALL Clot Types**






# Meet the growing NeVa Family

Introducing • **3 & 5 DROP ZONES:** multiple passes engineered into one  
• **LONGER WORKING LENGTHS:** extended thrombus interaction area

Vessel Diameters **2.0 - 3.5 mm**  
Ideal for MCA, ACA, PCA  
0.021 microcatheter



Vessel Diameters **2.0 - 4.5 mm**  
Ideal for ICA TIP, PROXIMAL MCA  
0.021 microcatheter



Vessel Diameters **3.5 - 6.0 mm**  
Ideal for PROXIMAL CAROTID ARTERY  
0.027 microcatheter

