


WANTY 035

OTW PTA Balloon Catheter



The BARTY MEDICAL PTA balloon catheters WANTY are indicated for peripheral transluminal angioplasty (PTA), for peripheral vascular system (including the iliac artery, femoral artery, popliteal artery)

Catheter Design	Coaxial
Balloon Design	Semi compliance
Nominal Pressure	6 ATM
RBP	Up to 16 ATM
Balloon Marker	2 swaged (zero profile)
Recommended Guidewire	0.035"
Sheath Compatibility	5F - 6F - 7F

MULTI SPECIFICATION AVAILABLE

Multi-size of balloon, diameter from 3.0 to 12.0 mm

LOW CROSS SECTION


Tip and folded balloon maintain a relatively low profile to get through and reduce damage to blood vessel

TRI-LAYER TUBING

Reducing push resistance of guidewire, providing support and pressure-bearing

MULTI FOLD

Fold based on balloon diameters, lowering profile to get through





WANTY 035

OTW PTA Balloon Catheter



PRODUCT SPECIFICATIONS

Shaft length (cm)	Diam (mm)	Balloon Length (mm)								
		20	40	60	80	100	120	150	200	220
80	3	H080030002035	H080030004035	H080030006035	H080030008035	H080030010035	H080030012035	H080030015035	-	-
	3.5	H080035002035	H080035004035	H080035006035	H080035008035	H080035010035	H080035012035	H080035015035	-	-
	4	H080040002035	H080040004035	H080040006035	H080040008035	H080040010035	H080040012035	H080040015035	H080040020035	H080040022035
	5	H080050002035	H080050004035	H080050006035	H080050008035	H080050010035	H080050012035	H080050015035	H080050020035	H080050022035
	6	H080060002035	H080060004035	H080060006035	H080060008035	H080060010035	H080060012035	H080060015035	H080060020035	H080060022035
	7	H080070002035	H080070004035	H080070006035	H080070008035	H080070010035	H080070012035	H080070015035	H080070020035	H080070022035
	8	H080080002035	H080080004035	H080080006035	H080080008035	-	-	-	-	-
	9	H080090002035	H080090004035	H080090006035	H080090008035	-	-	-	-	-
	10	H080100002035	H080100004035	-	-	-	-	-	-	-
12	H080120002035	H080120004035	-	-	-	-	-	-	-	
130	3	H130030002035	H130030004035	H130030006035	H130030008035	H130030010035	H130030012035	H130030015035	-	-
	3.5	H130035002035	H130035004035	H130035006035	H130035008035	H130035010035	H130035012035	H130035015035	-	-
	4	H130040002035	H130040004035	H130040006035	H130040008035	H130040010035	H130040012035	H130040015035	H130040020035	H130040022035
	5	H130050002035	H130050004035	H130050006035	H130050008035	H130050010035	H130050012035	H130050015035	H130050020035	H130050022035
	6	H130060002035	H130060004035	H130060006035	H130060008035	H130060010035	H130060012035	H130060015035	H130060020035	H130060022035
	7	H130070002035	H130070004035	H130070006035	H130070008035	H130070010035	H130070012035	H130070015035	H130070020035	H130070022035
	8	H130080002035	H130080004035	H130080006035	H130080008035	-	-	-	-	-
	9	H130090002035	H130090004035	H130090006035	H130090008035	-	-	-	-	-
	10	H130100002035	H130100004035	-	-	-	-	-	-	-
12	H130120002035	H130120004035	-	-	-	-	-	-	-	