

# iVascular<sup>®</sup>

therapies for living

# xperience

## Coronary balloon dilatation catheter

### Outstanding **deliverability**

Excellent pushability and flexibility to reach and cross the most challenging lesion

### Ultra-low tip and crossing profile

Optimized design

### Tungsten radiopaque markers

More flexibility with optimal trackability

## HYDRAX



Hydrophilic  
durable coating



[www.ivascular.global](http://www.ivascular.global)

## Real kissing balloon technique 6F compatible

Able to open the most challenging lesion.  
The low transition profile between the hypotube and the distal shaft allows the 6F kissing balloon technique



## Tungsten radiopaque markers

More flexibility with same visibility



# xperience features

- > **Rapid Exchange** catheter (RX)
- > Usable catheter length: **142 cm**
- > Compatible with **0.014"** guidewire
- > Compatible with **5F guiding catheter**
- > NP: **6 atm** | RBP: **16 atm** | ABP: **20 atm**
- > Tip profile: **0.016"**
- > Crossing profile: **from 0.021" up to 0.027"**
- > **2 polymeric radiopaque markers** with tungsten
- > **Semi-compliant**: 10-15%
- > **Deflation** time: 3s average



Product with CE mark, certified by Notified Body 0318

Balloon diameter (mm)	Usable catheter length: 142 cm					
	Balloon length (mm)					
	10	15	20	25	30	40
1.25	BC PRI4 150 125 010	BC PRI4 150 125 015	BC PRI4 150 125 020	-	-	-
	BC PRI4 150 125 110*	BC PRI4 150 125 115*	-	-	-	-
1.50	BC PRI4 150 150 010	BC PRI4 150 150 015	BC PRI4 150 150 020	-	BC PRI4 150 150 030	-
	BC PRI4 150 150 110*	BC PRI4 150 150 115*	-	-	-	-
2.00	BC PRI4 150 200 010	BC PRI4 150 200 015	BC PRI4 150 200 020	BC PRI4 150 200 025	BC PRI4 150 200 030	BC PRI4 150 200 040
2.25	BC PRI4 150 225 010	BC PRI4 150 225 015	BC PRI4 150 225 020	BC PRI4 150 225 025	BC PRI4 150 225 030	BC PRI4 150 225 040
2.50	BC PRI4 150 250 010	BC PRI4 150 250 015	BC PRI4 150 250 020	BC PRI4 150 250 025	BC PRI4 150 250 030	BC PRI4 150 250 040
2.75	BC PRI4 150 275 010	BC PRI4 150 275 015	BC PRI4 150 275 020	BC PRI4 150 275 025	BC PRI4 150 275 030	BC PRI4 150 275 040
3.00	BC PRI4 150 300 010	BC PRI4 150 300 015	BC PRI4 150 300 020	BC PRI4 150 300 025	BC PRI4 150 300 030	BC PRI4 150 300 040
3.25	BC PRI4 150 325 010	BC PRI4 150 325 015	BC PRI4 150 325 020	BC PRI4 150 325 025	BC PRI4 150 325 030	BC PRI4 150 325 040
3.50	BC PRI4 150 350 010	BC PRI4 150 350 015	BC PRI4 150 350 020	BC PRI4 150 350 025	BC PRI4 150 350 030	BC PRI4 150 350 040
3.75	BC PRI4 150 375 010	BC PRI4 150 375 015	BC PRI4 150 375 020	BC PRI4 150 375 025	BC PRI4 150 375 030	BC PRI4 150 375 040
4.00	BC PRI4 150 400 010	BC PRI4 150 400 015	BC PRI4 150 400 020	BC PRI4 150 400 025	BC PRI4 150 400 030	BC PRI4 150 400 040
4.50	BC PRI4 150 450 010	BC PRI4 150 450 015	BC PRI4 150 450 020	BC PRI4 150 450 025	BC PRI4 150 450 030	BC PRI4 150 450 040
5.00	BC PRI4 150 500 010	BC PRI4 150 500 015	BC PRI4 150 500 020	-	-	-

The availability of each reference for the sale is linked to the authorization of commercialization in the country of destination

\* Reference with one radiopaque marker centered in the balloon

Distributed by:



Manufactured by:

Life Vascular Devices Biotech S.L  
www.ivascular.global  
info@ivascular.global